

EasyDirect™ Moisture

Easy and Reliable Data Management v1.2



METTLER TOLEDO

Table of Contents

1	Introduction	3
1.1	End User License Agreement (EULA) and Disclaimer	3
2	Design and Function	4
2.1	EasyDirect™ Moisture Software	4
2.2	Functions	4
2.3	Installing EasyDirect™ Moisture.....	4
2.4	Configuring connection to EasyDirect™ Moisture in your Moisture Analyzer	6
3	Getting Started	8
3.1	General Settings	8
3.2	Check for Updates manually	8
3.3	Activating EasyDirect™ Moisture.....	8
4	Instrument Management	13
4.1	The Instrument View Window.....	13
4.2	Adding an Instrument.....	14
4.2.1	Connection via USB.....	15
4.2.2	Connection via RS232.....	15
4.2.3	Connection via Ethernet	16
4.2.4	Connection via Wi-Fi	16
4.3	Edit Instrument data and connection settings	17
4.4	Remove an Instrument	17
4.5	Date Filter.....	18
5	Test and Service Reminders	19
5.1	General Reminder Settings.....	19
5.2	Reminder settings for each instrument	20
5.3	Test and service status.....	21
6	Automatic Export	22
7	Measurement View	25
7.1	Transfer Measurement Results	25
7.2	Date Filter	27
7.3	Column Layout:	27
7.4	Measurement Results.....	27
8	Managing Results	28
8.1	Printing Results.....	28
8.2	Exporting Results	29
8.3	Deleting Measurement Results	29
8.4	Charts.....	30
8.4.1	Control Chart	30
8.4.2	Statistics Chart.....	31
9	Change Protection	32
10	Backup View	35
10.1	Import.....	35
10.2	Export.....	36
10.3	Edit Backup Properties	36
10.4	Delete Selected Backup	36
11	Data Transfer Indication	37

1 Introduction

This guide will help you to get started with EasyDirect™ Moisture v1.2 and explains the main features and functionalities of the software.

Easy and Reliable Data Management

Connect your HX204, HS153 or HC103 to EasyDirect Moisture software and improve your data management. Data collection starts at your PC, results are received even if the EasyDirect application is not open.

The EasyDirect Moisture software directly collects results and measured values of up to five instruments. Easily retrieve, manage, and filter results for each connected instrument.

Create charts based on results collected in EasyDirect, and print results on your local or network printer.

Your instrument results are stored and easily accessed using EasyDirect Moisture software. If you have any additional questions, METTLER TOLEDO is always available to assist you. Visit us on:

▶▶ www.mt.com/moisture

1.1 End User License Agreement (EULA) and Disclaimer

The software in this product is licensed under METTLER TOLEDO End User License Agreement (EULA) for software. When using this product you agree to the terms of the EULA.

▶▶ <https://www.mt.com/eula>

Disclaimer

This document may contain inaccuracies or errors of substantive or typographical nature. METTLER TOLEDO does not assume any responsibility or risk for your use of the information provided in this document. Nothing in this document shall be construed as a guarantee or warranty, either express or implied (including without limitation warranties as to merchantability or fitness for a particular purpose). The software described herein may be dependent on various external factors and configurations, such as system configuration, customer data, and operator control.

The implementation of the products this document relates to may vary consequently, the suitability of specific product configurations and applications must be determined by the customer and is not warranted by METTLER TOLEDO. If you rely on the information or any advice, opinion, or statement in this document, you are doing so at your sole risk. To allow improvements, the information in this document is subject to change.

2 Design and Function

2.1 EasyDirect™ Moisture Software

EasyDirect Moisture software consists of two parts:

- Background application
- Main user application

The background application runs automatically when the PC is running. The status will be indicated with an  icon in the Windows taskbar.

EasyDirect Moisture allows you to connect up to 5 moisture analyzers (HX204, HS153, HC103) via RS232 or USB interface and Ethernet or Wi-Fi .

The background application runs and logs results received from the connected instrument(s) regardless of whether you have EasyDirect Moisture open in a window or not. To view the results stored in the database, the main application has to be started.

2.2 Functions

The following represents a list of functions that can be performed using EasyDirect Moisture:

- Connecting the instrument through USB (Auto-Add)
- Manually adding the instrument through RS232
- Manually adding the instrument through Ethernet/Wi-Fi
- Viewing Instrument information
- Viewing measurement results
- Viewing tests and adjustment results
- Printing measurement results
- Exporting measurement results
- Automatic export of measurement & test results
- Remove instruments
- Change Protection
- View control charts and statistics charts
- Import Instrument backups, Method and settings files from USB drive
- Export Instrument backups, Method and settings files to USB drive
- Reminders for Tests and Service
- Automatic check for software updates

2.3 Installing EasyDirect™ Moisture


To install EasyDirect Moisture, a minimum of the PC requirements below should be satisfied:

PC requirements


Computer	Minimum	Optimum
RAM	4 GB	8 GB
Display Resolution	1024 × 768	1920 × 1080
Storage	5 GB free space	5 GB free space
Operating system	Minimum Windows 10 update 20H2	Windows 10 update 22H2

- 1 Download and run EasyDirect Moisture **Setup**.
 - ➔ The setup installs the Microsoft .Net framework version 4.8 as a pre-requisite before installing the application.
- 2 After successful installation of .Net Framework, restart the PC.
- 3 Select **Next** in the Welcome screen.

- 4 Read and accept the End-user License Agreement by clicking the corresponding checkbox.
- 5 Select **Next**.
 - ➔ The SQL server is upgraded based on the following scenarios:
 - a. EasyDirect Moisture first-time installation: the setup checks whether there is an instance already of name EasyDirect with SQL server version 2022. If there is none, then SQL server 2022 is installed.
 - b. Lower version of EasyDirect Moisture is already installed:
 - If the existing SQL server version is 2019, the application will automatically upgrade the SQL server from 2019 to 2022 along with the application upgrade.
 - If the existing SQL server version is not compatible, the following message appears and installation fails. You must upgrade manually the SQL server 2012 from SP1 to SP4 and restart the PC. Run the setup again to upgrade the SQL server from 2012 to 2022 along with the application upgrade.



- 6 Optionally, select the installation path.
- 7 Click **Install** button.
- 8 The setup proceeds with the product installation.
- ➔ On successful installation, the screen below is displayed.



➔ After you restart your PC, the Welcome screen is displayed.

2.4 Configuring connection to EasyDirect™ Moisture in your Moisture Analyzer

Instrument requirements

Moisture Analyzer	Instrument software
HX204	v3.31 (56358) and v2.50 (53385)
HS153	v3.31 (56358) and v2.50 (53385)
HC103	v2.11 (56358) and v1.70 (53385)

To configure the USB connection:

- 1 Power on your Moisture Analyzer [HX204/HS153/HC103].
- 2 From the Home screen, select **Settings**.
- 3 Select **Instrument Settings**.
- 4 Select **Peripheral**.
- 5 Select **New...** to add a peripheral.
- 6 Select EasyDirect Moisture.
- 7 To connect EasyDirect Moisture instruments, Choose **USB** as device port for connection using a USB cable.
- or -
Choose **RS232C** as device port for connection using an RS232 cable.
- 8 Save the peripheral.
- 9 Enable the peripheral.

To configure the Ethernet or Wi-Fi connection:

- 1 Power on your Moisture Analyzer [HX204/HS153].
- 2 From the Home screen, select **Settings**.
- 3 Select **Instrument Settings**.
- 4 Select **Peripheral**.
- 5 Select **New...** to add a peripheral.
- 6 Select EasyDirect Moisture.
- 7 Select device interface as **Network** and default port is **8100**.
- 8 Save the peripheral.
- 9 Enable the peripheral.

Note


Configure Ethernet or Wi-Fi (WLAN) through **Network interfaces** menu under **Instrument settings**.

3 Getting Started

- In the Windows **Start** menu, select METTLER TOLEDO > EasyDirect Moisture
 - OR -
 - double-click on the  shortcut on the desktop.
- ➔ EasyDirect Moisture is started.

3.1 General Settings

- 1 Select  **Options** from the ribbon.
 - ➔ The Options dialog appears:




- 2 Select the desired **Language** from the drop-down list. The default language is set to English. EasyDirect Moisture supports the following languages: English, German, French, Italian, Spanish, Russian, Polish, Czech, Hungarian, Portuguese, Chinese, Japanese, and Korean.
- 3 By unchecking **Show notifications in tray**, you will not receive tray notifications from EasyDirect Moisture.
- 4 If you do not want EasyDirect Moisture to check for updates automatically, you can uncheck **Check for updates automatically**.
- 5 **Save** the settings to reflect the changes.

3.2 Check for Updates manually

The application automatically checks for new updates every day when this option is selected in **General settings** and downloads the latest version of the application. You can also manually check for updates from the  **About** window.

EasyDirect Moisture Update Intimation

EasyDirect Moisture provides a new software update in the About view. When you click About, a new window appears and you can click check for updates.




New update on EasyDirect Moisture is found during the following conditions:

- 1 When user manually tries to check for a new updates.
- 2 During EasyDirect Moisture launch.

3.3 Activating EasyDirect™ Moisture


The trial version of EasyDirect Moisture can be used for 45 days without a license. After this, a license is required for continued use of the software. A license can be purchased through your METTLER TOLEDO sales-person.

▶▶ <http://www.mt.com/contacts>



Steps to activate EasyDirect™ Moisture product:


- 1 Start your installed EasyDirect™ Moisture application and note / copy the License Registration Information (LRI).




- 2 Sign in to <https://activation.mt.com> and refer to **Activation**.
- 3 Enter the EasyDirect™ Moisture license key and enter / paste the LRI and click submit.




Activation / Reactivation / System Transfer



- The Software Portal provides a Short Activation Code (SAC) in a pop-up box. You can copy the SAC displayed on the pop-up window.




- Enter or paste the SAC into your EasyDirect™ Moisture activation window and click **Register**.



- ➔ EasyDirect™ Moisture application is now successfully activated.

View activated system on Software Portal



To add or edit system information

- 1 Select an activated system and click the magnifier.

Activation / Reactivation / System Transfer

Product Subname Username

In order to activate your product, first export the license key file using the "License manager" in your product. The license key file contains the hardware ID of your PC (Server) and all the license keys of your installation. On this page, browse for the path to the exported license key file (.lic) and click on "Submit" to upload.

License key file: No file chosen

Add license key + License Registration Information (LRI)

Add a single license key of a simple SW product or add a Transfer key to receive an activated system from another user account.

License key: LRI:

Add Transfer Key

Add a transfer key to receive an activated system from another user account.

Transfer Key:

Activated products Non-Activated products

Please select to change, update or view details of a system.

Product Subname	System name	Server name	Server location	Hardware ID / LRI
EasyDirect Moisture				LEBZ-2595

- 2 Select **Edit System** in License overview and add/change system information.

METTLER TOLEDO

ACTIVATION PORTAL

[Home](#) [Knowledge Base](#) [Activation](#) [Profile](#) [Register](#)

License overview

Product Subname EasyDirect Moisture

System Name [Set Transfer Key](#)

Server Name [Edit System](#)

Server Location

Hardware id / LRI LEBZ-2595

Activated editions

You can reactivate your license, by clicking on "Reactivate" button.

Product subname	License Key	Uploaded	Activated	Short Activation Code
EasyDirect Moisture	6R6L6-K4NJ7-65622-62PJS-6UQPR	6/28/2022	6/28/2022	LE3D2-S7J8L-66JEV

- 3 You can view the activated systems including system information in the activation portal.

Activation / Reactivation / System Transfer

Product Subname

Username

In order to activate your product, first export the license key file using the "License manager" in your product. The license key file contains the hardware ID of your PC (Server) and all the license keys of your installation. On this page, browse for the path to the exported license key file (.lic) and click on "Submit" to upload.

License key file: No file chosen

Add license key + License Registration Information (LRI)

Add a single license key of a simple SW product or add a Transfer key to receive an activated system from another user account.

License key: LRI:

Add Transfer Key

Add a transfer key to receive an activated system from another user account.

Transfer Key:

Activated products

Non-Activated products

Please select  to change, update or view details of a system.

Product Subname	System name	Server name	Server location	Hardware ID / LRI
EasyDirect Moisture				1EBZ-2595


4 Instrument Management

4.1 The Instrument View Window

Navigation:  **Instrument View**

For instrument requirements see [Configuring connection to EasyDirect™ Moisture in your Moisture Analyzer ▶ Page 6].


The instrument view provides an overview of instrument details such as test and service status.



Connection Status indication

	Green	The instrument is successfully connected and is monitored.
	Grey	The instrument is added but is currently disconnected
	Red	Error in connection

Expand an instrument by clicking the  icon or double click to view the instrument details



Expanded view displays details about the instrument including, Terminal Software, Drying Unit Software, Hardware details, Connection settings, Company Information and Test and Service status.

The **Test and Adjustment results (Weight, Temperature and SmartCal)** can be viewed here.

Connection Type Indication

EasyDirect Moisture supports every user to identify the different connection types along the instrument results in the **Instrument view**.


	The instrument is successfully connected via USB.
	The instrument is successfully connected via RS232.
	The instrument is successfully connected via Ethernet/Wi-Fi.

Note

In case of application restart, connection type indication of the disconnected USB instrument will be changed to RS232 icon. It reflects USB icon once the instrument is connected again through USB.

4.2 Adding an Instrument

The software can handle up to five connected instruments simultaneously. The software detects connection of a sixth instrument, but no monitoring or collection of the new instrument data will occur. This attempt to connect a sixth instrument also results in an error message:



- To add a new instrument, you have to remove an existing instrument from the list.


When the instrument is added to EasyDirect Moisture, you can disconnect and re-connect from EasyDirect Moisture using  **Connect** or  **Disconnect** in the ribbon.

4.2.1 Connection via USB


When an instrument is connected to the PC via USB, it is automatically detected. If the instrument is found to be compatible, a message box is displayed, and you can add it to the list of monitored instruments. Once added to this list, the application will automatically connect to the instrument when the instrument is re-connected and switched on (unless the instrument is subsequently removed from the list).

For Moisture Analyzers connected via USB cable, use the following procedure:

- 1 Switch the Moisture Analyzers on.
- 2 Connect a Moisture Analyzer to the PC using a USB cable.
 - ➔ EasyDirect Moisture should start automatically. The message box below is displayed:



- 3 Click **Yes** to add the detected Moisture Analyzer.
 - ➔ Once the Moisture Analyzer is added, you can modify the auto generated System ID for the instrument in the Instrument data popup.



- ➔ Any measurement or test performed in the Moisture Analyzer after adding will be transferred automatically to EasyDirect Moisture


4.2.2 Connection via RS232

For the instruments connected through the RS232 port, you have to manually add the instrument by selecting the correct COM port for the connection.

If you are using the RS232 port to connect to the Moisture Analyzer, use the **Add Instrument** window to add the instrument manually:

Navigation: Manage

- 1 Select **Add instrument**.



- 2 In the **Add Instrument** window, choose the instrument type and choose the RS232 COM port.

4.2.3 Connection via Ethernet

Ethernet connectivity provides user to connect the instruments in easier way from working environment. This mode of connectivity is supported only for **HX204** and **HS153** type instruments with the firmware versions v3.31 (56358) and v2.50 (53385). By default, the device port field is pre-filled with the value 8100.



- 1 Click on **Manage menu**, select **Add instrument**.
- 2 In Add Instrument window, choose the interface type as **Ethernet**.
- 3 You can manually enter the IP address and Device Port number.
- 4 Click **Connect** to start the connection through Ethernet connectivity.

4.2.4 Connection via Wi-Fi

For the instruments connected through Wi-Fi, you are allowed to connect only by using **LM007 Wi-Fi dongle**.



- 1 Click on **Manage menu**, select **Add Instrument**.
- 2 In **Add Instrument** window, choose the interface type as **Wi-Fi**.

- 3 You can manually enter the IP address and Device Port number.
- 4 Click **Connect** to start the connection through Wi-Fi connectivity.

Note

- To know more information about LM007 Wi-Fi dongle, you can contact METTLER TOLEDO sales team.

4.3 Edit Instrument data and connection settings

Navigation:  **Manage**

You can edit the auto generated System ID, Company name, Department, Location and the user who is responsible for the instrument by selecting **Edit selected instrument data** in the **Manage** menu.


To change the RS232 COM port and connection settings, select **Edit selected instrument connection settings** in the **Manage** menu.

4.4 Remove an Instrument

Navigation:  **Manage**


- 1 Select the instrument to be removed from the application.
- 2 Disconnect the instrument by selecting  **Disconnect** from the ribbon.
- 3 Select  **Manage** from the ribbon.
- 4 Click **Remove instrument**.

➔ The following message appears:



- 5 Click **Yes**.

➔ The following message appears:




- 6 Click **Yes** and choose a folder to where the measurement results, test results and instrument backups will be exported. Once exported the instrument will be removed.
-or-
Click **No** to remove the selected instrument without exporting data.
- or -
Click **Cancel** to cancel the remove operation. In this case, the selected instrument will not be removed.

4.5 Date Filter

Navigation:  **Filter**

To filter Tests and adjustments by date range, select the **filter option**  in the ribbon.



The following Date filters can be applied on the **Tests and Adjustments** results displayed in the expanded instrument view:

- **No date filter** - displays all the Test and Adjustment results without filtering
- **Date range** - Filter results based on Date range
- **Last No. of days** – Filter results based on the last number of days
- **Today** - Filter and display only today's Tests and Adjustments


5 Test and Service Reminders

Navigation:  Reminders

Test and Service Reminders allow you to define the frequency for routine testing of balance (weight), heating unit (temperature) and overall instrument test (SmartCal) and benefit from the integrated reminder function. You can schedule reminders for tests and service that have to be performed on your Moisture Analyzers. Default reminders can be set to all connected instruments using **General Reminder Settings** or to specific instruments using **Selected instrument reminders**.

5.1 General Reminder Settings

General reminder settings can be accessed from the  **Reminders** menu.



General reminder settings

Weight test interval
Define weight test interval: 30 day(s)
 Activate status and reminder


Temperature test interval
Define temperature test interval: 30 day(s)
 Activate status and reminder

SmartCal test interval
Define SmartCal test interval: 30 day(s)
 Activate status and reminder

Service reminder
Define service reminder: 30 day(s)
 Activate status and reminder

Save Cancel

Default reminders can be set for the Weight test, Temperature test, SmartCal test and Service of all connected instruments. You can select the test and service interval. You can deactivate the reminder by unchecking **Activate Status and reminder**. When the tests or Service of the connected instruments are due, you will be notified with a reminder message.




Note

- On saving the **General reminder settings**, the reminders for all the instruments will be overwritten.
- The Service date is received from the instrument (it is set for the first time in production and then a MT service engineer can re-define it with LARS)

5.2 Reminder settings for each instrument

Reminder settings of an instrument can be changed by selecting the **Selected instrument reminder** under the  **Reminders** menu.



Weight test interval

Last weight test date: 2/21/2019

Define weight test interval: 60 day(s)

Next weight test date: 4/22/2019

Activate status and reminder

Temperature test interval

Last temperature test date: 2/21/2019

Define temperature test interval: 60 day(s)

Next temperature test date: 4/22/2019

Activate status and reminder

SmartCal test interval

Last SmartCal test date: 3/29/2019

Define SmartCal test interval: 14 day(s)

Next SmartCal test date: 4/12/2019

Activate status and reminder

Service reminder

Next service date: 3/30/2019

Define service reminder: 30 day(s)

Activate status and reminder

Save Cancel

You can define the Weight test, Temperature test, SmartCal test. The service reminder defines the days you will be reminded before this service due date. When you perform a test in the instrument, the results are automatically received by EasyDirect Moisture. Based on the last test date and the defined interval, the next test dates are computed. The service and test status are also displayed in the expanded view of the instrument in the instrument view. You can uncheck Activate Status and reminder, if you do not wish to see the reminders and status.

Note

- If you have not performed any test yet, the next test date can be configured manually.
- The last test date and next test dates are updated when a new test result is received from the Moisture Analyzer.

5.3 Test and service status

Test and service status of the connected instruments can be observed in the expanded view of the instrument.

The screenshot displays the METTLER TOLEDO EasyDirect™ Moisture software interface. The top menu bar includes options like Instrument view, Measurement view, Backup view, Disconnect, Manage, Automatic export, Reminders, Filter, Export results, Print results, Options, Help, and About. Below the menu is a table of instruments with columns for Status, System ID, Type, Location, Instrument status, Weight test status, Temp. test status, SmartCal test status, and Service status. Two instruments are listed: HC103 1 (HC103, Ready) and HX204 1 (HX204, Ready). The HX204 1 instrument is expanded to show detailed information:

- Instrument:** System ID: HX204 1, Type: HX204
- Terminal software:** Identification: 30535417Dtest, Version: 3.30.01, TDNR: 32.28.0.389.142
- Drying unit software:** Identification: 11671197L, Version: 1.60D, TDNR: 51.28.5.3617.1195
- Hardware:** Terminal SNR: 9876543210, Drying unit SNR: 9876543210
- Company information:** Company name, Department, Location, Responsible
- Connection settings:** IP address: 127.0.0.1, Device port: 8100
- Test and Service status:**
 - Weight test status: Test overdue for 19 day(s).
 - Temp. test status: Test interval not configured.
 - SmartCal test status: Test interval not configured.
 - Service status: Next service date not received.

Below the instrument details is a 'Tests and adjustments' section with tabs for Weight, Temperature, and SmartCal. A table shows test history:

Date/Time	System ID	Mode of test	Status	User	Nominal weight	Actual Weight
08-Dec-20 12:17:42 AM	HX204 1	Weight test internal				
16-Nov-20 07:36:00 PM	HX204 1	Weight adjustment...	Done	System	100.000 %	
16-Nov-20 07:03:00 PM	HX204 1	Weight adjustment...	Done	System	100.000 %	
20-Oct-20 12:17:00 PM	HX204 1	Weight test external	Passed	System	100.000 g	100.001 g

Selected row count: 1. Auto export status:

Test status indication

	Test reminder not configured (or) Activate Status and Reminder unchecked in Reminder settings (or) Last test not received
	Last Test Passed, ready for operation
	Test due today
	Last Test Failed (or) Test overdue for x days

Service status indication

	Activate Status and Reminder unchecked in Reminder settings (or) Next service date not received
	Service done
	Service due in x days (or) Service due today
	Service overdue for x days

6 Automatic Export

- 1 To schedule an automatic export across all the connected instruments, select **Global auto export**. For any instrument specific export, choose **Selected instrument auto export**. Automatic export supports only for **HX204** and **HS153** type instruments with the firmware versions v3.31 (56358) and v2.50 (53385).

Global auto export

Export On

Export by

Immediate

Every 4 Hour Start time 20/12/21 12:51:40 PM

Last configured start time 20/12/18 10:37:02 AM

Last exported time 20/12/21 10:37:02 AM

Measurements and tests Measurement;Weight;Temper...

Export file type PDF

File name prefix AE

Target folder C:\Users\Public\METTLER TOLE ...

One result in a file


All results of a category in a file

Please note that the maximum number of records for PDF format, that can be exported to a single file is 2000. After that the records will be exported to another file.

Save Cancel

- 2 Switch **Export ON**.
- 3 If already configured, last configured start time and last exported time will be displayed.
- 4 If you want to export as soon as the data is received, choose **Immediate**.
- 5 If the export should happen in certain frequency, choose **Every – hour** and set the hour frequency and start time. If you want to export once a day, please set the hour frequency to 24 hours.
- 6 Here auto export will start exporting the data only from the selected start time and does not consider other results.

7 You can choose the measurements including test and adjustments to export.




8 Export can be done in CSV, XML and PDF formats.

9 Each measurement will be exported to separate file if 'One result in a file' option is chosen.

10 All results of each category will be exported to a single file, if the **results of each category in a file** is chosen.


11 If you select **CSV** as the file type, then by default **One result in a file** option is selected. Here in CSV file type, every results are stored in a separate file.



12 Choose the folder for the exported file to be saved in **Target folder**. Ensure to choose public folders to avoid access errors.

Auto export Status

Once you switched into auto export feature, you will find the present status on the bottom of your screen. Every status is differentiated through Grey, Green and Red color.



	Grey	Auto export is disabled.
	Green	Enabled for any of the instruments and successfully exported.
	Red	Auto export is unsuccessful for any of the instruments.


Note

- The maximum number of records for PDF format that can be exported to a single file is 2000. After that the records will be exported to another file.
- You will find the date and time of a last auto export in the mouse too tip, irrespective of a successful or unsuccessful completion.
- Auto export status is displayed in generic and not displayed as instrument specific.
- In case of any connectivity issue, you can manually transfer the result by using **Export** or **Print** button.

7 Measurement View

Navigation: Measurement view

The Measurement view lists the measurement results received from all the connected instruments. You can select multiple results at a time. To do this, hold the Shift key and click on results to select a range of results, or hold the Ctrl key and click on the individual results you wish to select. This is useful for printing or exporting results.



Date/Time	System ID	User name	Method name	End result	Total
01.04.2019 16:40:00	HX204 2	Administrator	POLENTA QP	12.43 %MC	15
29.03.2019 15:48:00	HS153 - 01	Administrator	MILKPOWDER	3.79 %MC	0:5
21.03.2019 16:54:00	HX204 1	Administrator	TEST ID	4.57 %MC	3:1
12.03.2019 09:10:00	HX204 2	Administrator	HAZELNUT	4.56 %MC	6:1
25.02.2019 15:58:00	HS153 - 01	Administrator	MILKPOWDER	6.17 %MC	3:0
21.02.2019 17:11:00	HS153 - 01	Administrator	MILKPOWDER	5.13 %MC	8:4
21.02.2019 16:47:00	HS153 - 01	Administrator	MILKPOWDER	2.97 %MC	7:2
21.02.2019 16:36:00	HS153 - 01	Administrator	MILKPOWDER	5.91 %MC	8:1

Duration	Result value	Result unit
00:00:00	0.00	%MC
00:00:10	0.11	%MC
00:00:20	0.70	%MC
00:00:30	1.66	%MC
00:00:40	2.76	%MC
00:00:50	3.85	%MC
00:01:00	4.80	%MC
00:01:10	5.53	%MC
00:01:20	6.10	%MC
00:01:30	6.59	%MC
00:01:40	6.97	%MC
00:01:50	7.35	%MC

7.1 Transfer Measurement Results

Transferring results automatically to EasyDirect Moisture

Perform an analysis using the connected moisture analyzer. Wait for the results of measurement to be transferred to EasyDirect Moisture. Measurement results will be transferred automatically from the analyzer to EasyDirect Moisture.

Measurement results from the instrument can be viewed in the **Measurement** view.

Transferring results manually to EasyDirect Moisture

When Instrument is not connected with EasyDirect Moisture, results are manually transferred through following procedure.

- 1 To transfer the results, Go to **Results** and click Method name in the Instruments.

2 **Click** Performed result and **Select** Print image button.

Print result	
Printer	EasyDirect Moisture
Printout type	Standard >
Printout language	English >
Include intermediate results	OFF >

Cancel Print

3 **Click** Print button and your results will be transferred to EasyDirect Moisture.

Measurement results from the instrument can be viewed in the **Measurement** view.

Measurement view offers the following features:

- Search for a specific text in the list - If you enter a search term in **Search box**, the rows will be filtered to those containing the search term.
- **Resize** the columns
- **Re-Order** the columns by dragging them
- **Sort Date and Time** by clicking the header
- **Filter** results by clicking filter icon in the column header
- Choose one of the pre-defined column layout or save your own layout

7.2 Date Filter

To filter the measurement results by date range, select  **Filter option** in the ribbon.

The following Date filters can be applied on the measurement result:

- **No date filter** - displays all the Test and Adjustment results without filtering
- **Date range** - Filter results based on Date range
- **Last No. of days** – Filter results based on the last number of days
- **Today** - Filter and display only today's Tests and Adjustments

7.3 Column Layout:


EasyDirect Moisture provides the following predefined column layouts for the measurement view:

- Default layout
- All Columns – Displays all columns available

You can create a customized layout by resizing, reordering or adding/removing columns to the layout and save it as new layout. You can choose any of the available layouts in the **Column layout** drop-down in the measurement view.

You can re-order the columns by dragging the column headers. To add/remove columns to the layout, right click on the header row and choose **Column Chooser**. The **Column Chooser** lists all the available columns. You can include a column by checking the column.

You can save the customized layout by selecting **Save as new layout** in  **Column layout**.



Column Layout provides the options below:

- **Save current layout** – After customizing a layout, you can save it. A predefined layout cannot be saved. After making any changes to the predefined layout, you can save it as a new layout.
- **Delete current layout** – Delete the selected layout. A predefined layout cannot be deleted.

7.4 Measurement Results

The measurement list displays results based on the date filter. When you select a result, the **Drying curve** will be displayed to the right. **Intermediate values** from the measurement will be displayed below the curve.


8 Managing Results

8.1 Printing Results


Navigation: Print results

Measurement results from Measurement view and Test and Adjustment results from Instrument view can be printed to a connected printer using the **Print results** option. Up to 150 results can be printed at once by EasyDirect Moisture.

- 1 Select the results you want to print.
- 2 Click  **Print results**.
 - ➔ **The Print results dialog opens.**



- 3 Select any one of the options.
- 4 Click **Print**.
 - ➔ A print preview is displayed




- 5 Select **Print** to get the results printed.

8.2 Exporting Results

Navigation: **Export results**

Measurement results from Measurement view and Test and Adjustment results from Instrument view can be exported using the **Export results** option. You can export the results as CSV, PDF or XML format.

- 1 Select the results you want to export.
- 2 Click  **Export results**.
 - ➔ The **Export results** dialog opens.



- 3 Select the **Export file type** as CSV, PDF or XML.
- 4 Enter the **File name prefix**, if you want a prefix for the exported file names.
- 5 Select the destination folder in **Export folder**.
- 6 Click **Export**.
 - ➔ The results are exported and saved as .csv, .pdf or .xml files (based on the file type selected) in the specified folder. Each result is exported as a separate file.

8.3 Deleting Measurement Results

Navigation: **Delete results**

Only measurement results can be deleted, while the Test and Adjustment results in the Instrument view cannot be deleted.

- 1 Select the rows you want to delete from the Measurement view.
- 2 Click  **Delete results** from measurement view window.
 - ➔ The selected measurement results are deleted after confirmation.

8.4 Charts

Navigation: Charts

Measurement results from Measurement view and SmartCal results of an instrument from Instrument view can be visualized graphically. Two types of charts can be plotted with the results in EasyDirect Moisture:


- Control Charts
- Statistics Charts

8.4.1 Control Chart

The **Control Chart** is used to visualize the results over time with control limits defined in the method.

You can plot Control charts by selecting the **Control chart** option from the **Charts** menu.

- 1 Filter the Results for same method name by clicking on the filter icon next to the **Method name** column header and choosing only the method you want to plot. Ensure all the results have the same unit.
 - 2 Results can be excluded from being plotted by unchecking them.
- ➔ Test Measurement results and Aborted measurement results will be excluded in the chart.




The control chart displays the control limits based on the limits set for the most recent selected measurement. If the method definition is modified, the results that belong to the latest method definition are plotted with solid marker point and the results for the older method definition are plotted with outlined marker points.

You can use the **Copy chart to clipboard** button to copy the chart to the clipboard or use the **Save chart as image** button to save it as an image file.

8.4.2 Statistics Chart

The **Statistics Chart** is used to visualize the results over time including visualization of mean and standard deviation of the selected results. You can plot Statistics charts by selecting the **Statistics chart** option from the **Charts** menu

- 1 Filter the Results for same method name by clicking on the filter icon next to the **Method name** column header and choosing only the method you want to plot. Ensure all the results have the same unit.
 - 2 Results can be excluded from being plotted by unchecking them.
 - ➔ Test Measurement results and Aborted measurement results will be excluded in the chart.
- ➔ Mean and the Standard deviation of the selected results are displayed:



9 Change Protection

Navigation: Options

Change protection feature enables you to configure permissions for the protected mode. This feature provides you options to allow or deny access to:


- Delete measurement results
- Add or remove instruments
- Backup/Restore instrument data
- Edit/Modify the scheduler etc.

Initially, **Change protection** will be disabled and you can access all the features. You have to enable Change protection to limit access in the protected mode.

To enable protection, follow the steps below:

- 1 Select the **Set Protection Password** option from the **Options** menu.

➔ Set Protection Password window will appear




The 'Set protection password' dialog box contains three input fields: 'Enter password' with four dots, 'Re-Enter password' with four dots, and 'Password recovery e-mail address' with the text 'user@company.com'. At the bottom are 'Save' and 'Cancel' buttons.

- 2 Enter a password and confirm the password.
- 3 Enter a valid e-mail address to recover the password later if forgotten in the **Password recovery e-mail address** field.


- 4 **Save** the password.

➔ A window prompting to configure Mail server Settings will appear.



The 'EasyDirect™ Moisture' dialog box features a green question mark icon and the text 'Do you want to configure Mail server settings now?'. It has 'Yes' and 'No' buttons at the bottom.

- 5 Click on **Yes** to configure **Mail server Settings**
- 6 Click No to configure later by selecting Mail server Settings from the options menu.



The 'E-mail server settings' dialog box has two sections. The top section, 'E-mail server settings', includes: 'SMTP Host *' (outlook.mt.com), 'Port *' (25), 'Time out (seconds) *' (20), 'Use SSL' (checkbox), 'Authentication type *' (Windows dropdown), and 'Sender e-mail address *' (support@company.com). The bottom section, 'E-mail settings verification', includes: 'Test e-mail address' (user@company.com) and a 'Test e-mail' button. 'Save' and 'Cancel' buttons are at the bottom.

- 7 Fill in the **E-mail server settings** window. If you are not sure what to fill in, please contact your local IT.

- 8 To test the Mail server settings, enter a valid e-mail id in the **Test e-mail address** field and use the **Test e-mail** button to verify if a test e-mail is received.


9 **Save** the settings.

➔ Now, you have enabled protection.

To limit the features you have to enter the unprotected mode.

1 Select the **Enter Unprotected Mode** from the **Options** menu.

➔ Enter Unprotected Mode window will appear.



2 Enter the protection password in the **Password** field and **Log in**.

➔ Application enters Unprotected Mode on successful login.


3 In case you forgot the password, use the **Forgot Password** link.

➔ The protection password will be sent to the **Password recovery e-mail address** set while enabling protection.

➔ After successful login, you are in unprotected mode. Now, you can define permissions for the protected mode.

4 Select the **Protection Settings** from the **Options** menu.

➔ The Protection settings window will appear.



You can choose the permissions for various features by checking them. Unchecked features will not be accessible in the protected mode.

Change password enables to change the protection password

You can also change the **Password recovery e-mail address**

Clicking **Remove Protection** will remove the password and protection settings and you will have unlimited access to all the features. Mail server settings will not be removed on removing protection.


After configuring the permissions, you can return to the protected mode.

- Select **Leave unprotected mode** from Options menu.
- ➔ Permissions selected in the protection settings will be applied.

10 Backup View

Navigation: Backup View


An instrument backup is vital in case of instrument failure. Use EasyDirect Moisture to store backups and instrument settings. **Backup view** enables you to import the instrument backups, method and settings files from a USB drive to EasyDirect Moisture and export to the USB drive.



10.1 Import

To import a backup,

- 1 Choose the **Instrument backups** or **Method and settings files** based on what you want to import.
- 2 Select  **Import** from the ribbon in the **Backup View** screen.
 - ➔ The Import backup window will be displayed.



- 3 Browse for the backup file or the Methods and settings file in the USB drive.

- 4 Click **Save**.
- ➔ Once saved, the Backup details will be displayed in the Backup view.

Note

- You can import a backup of an instrument only if the instrument is added to EasyDirect Moisture.

10.2 Export

- 1 Select any backup from **Instrument backups** or the **Method and settings files** tab in the Backup view.
- 2 Select **Export** from the ribbon in the Backup View.
 - ➔ The Export backup window will be displayed.



- 3 Browse for the file in the USB Drive to restore the backup.
- 4 Click **Restore**.
 - ➔ Backup will be stored in the selected folder successfully.

10.3 Edit Backup Properties

Backup Properties can be edited using the **Edit Backup Properties** option. This can be accessed from the **Manage** menu in the Backup view screen.

- 1 Select any backup from the **Instrument backups** or **Methods and Settings files** tab in the Backup View screen.
- 2 Choose the **Edit Backup Properties** option from the **Manage** option.
 - ➔ The Edit Backup Properties window will appear.
- 3 Edit the Backup name or the Comments of the Backup.
- 4 Click **Save**.

10.4 Delete Selected Backup

A Backup can be deleted using the **Delete Selected Backup** option.

- 1 Select any backup.
- 2 Choose the **Delete Selected Backup** option from the **Manage** menu to delete the backup.

11 Data Transfer Indication

User finds a visual indication in the screen when any data is transferred from the instruments to EasyDirect Moisture.

The screenshot displays the METTLER TOLEDO EasyDirect™ Moisture application interface. The top menu bar includes options like Instrument view, Measurement view, Backup view, Disconnect, Manage, Automatic export, Reminders, Filter, Export results, Print results, Options, Help, and About. Below the menu is a table of instruments with columns for Status, System ID, Type, Location, Instrument status, Weight test status, Temp. test status, SmartCal test status, and Service status. The selected instrument is HX204 1, which is in a 'Ready' status. Below the table, there are sections for Instrument, Company information, Terminal software, Drying unit software, and Hardware. The 'Tests and adjustments' section is active, showing a table of test results. At the bottom, a progress bar indicates 'Transfer in progress' with a green bar and a 'Selected row count: 1' label.

Status	System ID	Type	Location	Instrument status	Weight test status	Temp. test status	SmartCal test status	Service status
Ready	HC103 1	HC103		Ready				
Ready	HX204 1	HX204		Ready				

Date/Time	System ID	Mode of test	Status	User	Nominal weight	Actual Weight
08-Dec-20 12:17:42 AM	HX204 1	Weight test internal				
16-Nov-20 07:36:00 PM	HX204 1	Weight adjustment...	Done	System	100.000 %	
16-Nov-20 07:03:00 PM	HX204 1	Weight adjustment...	Done	System	100.000 %	
20-Oct-20 12:17:00 PM	HX204 1	Weight test external	Passed	System	100.000 g	100.001 g

- 1 When there is any data transfer from the instruments to EasyDirect Moisture application, green indication blinks on the screen.
- 2 Additionally, you will find a progress bar with the message **Transfer in progress** with a green indication. This message applies to all the data transfer irrespective of views.
- 3 When you transfer a data which is already available in EasyDirect Moisture data list, a progress bar appears on the screen but the duplicate data is not added to the list.

To protect your product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of this product for years to come.

Please request full details about our attractive terms of service.

► www.mt.com/service

www.mt.com

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contacts

Subject to technical changes.

© 04/2023 METTLER TOLEDO. All rights reserved.