WMF

English	User Manual WMF Weigh Modules
Deutsch	Benutzerhandbuch WMF Wägemodule
Français	Guide de l'utilisateur WMF Modules de pesage
Español	Manual de usuario WMF Módulos de pesaje
Italiano	Manuale per l'utente WMF Moduli di pesata
Nederlands	Handleiding WMF Weegmodules
Português	Manual do Usuário WMF Módulos de Pesagem

METTLER TOLEDO

User Manual **WMF**

English

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Deutsch

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Français

Manual de usuario **WMF**

Español

Manuale per l'utente **WMF**

Italiano

Handleiding **WMF**

Nederlands

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Português

Table of Contents

1	Introduction	3
1.1	This User Manual.....	3
1.2	Further Documents.....	3
1.3	Manufacturer Information.....	3
2	Safety Information	4
2.1	Intended Use.....	4
2.2	Definition of Signal Words and Warning Symbols.....	4
2.3	Product Specific Safety Notes.....	4
3	WMF Weigh Modules	5
3.1	Specifications.....	5
3.2	Overview.....	6
3.3	Scope of Delivery.....	7
4	Operation	8
4.1	Electrical Connections.....	8
4.2	Module Status LED.....	8
4.3	Environmental Conditions.....	8
4.4	Ingress Protection.....	8
4.5	Commissioning.....	9
5	Maintenance	10
5.1	Cleaning.....	10
5.2	Calibration and Adjustment.....	10
5.3	Disposal.....	10

1 Introduction

1.1 This User Manual

This User manual contains all information for the **operator** of the product.

- Read this User manual carefully before use.
- Keep this User manual for future reference.
- Pass this User manual to any further owner or user of the product.

1.2 Further Documents

In addition to this printed User manual you can download the following documents from Documentation WMF

► <http://www.mt.com/WMF>

- Installation Manual
- SAI Reference Manual
- MT-SICS Reference Manual
- Technical data sheet

1.3 Manufacturer Information

The contact information of the manufacturer of the product is as follows:

- **Name:** METTLER-TOLEDO GmbH
- **Weblink:** <http://www.mt.com>
- **Physical address:** Im Langacher 44, 8606 Greifensee, Switzerland

2 Safety Information

2.1 Intended Use

- Use the product only for weighing in accordance with this User manual.
- The weigh module is intended for indoor use only.
- Any other type of use and operation beyond the limits of technical specifications is considered as not intended.

2.2 Definition of Signal Words and Warning Symbols

Safety notes are marked with signal words and warning symbols. These show safety issues and warnings. Ignoring the safety notes may lead to personal injury, damage to the load cell, malfunctions and false results.

Signal Words

CAUTION	Hazardous situation with low risk, resulting in damage to the device or the property or in loss of data, or minor or medium injuries if not avoided.
Attention	Important information about the product (no symbol)
Note	Useful information about the product (no symbol)

Warning Symbols

General hazard

Electrical shock

2.3 Product Specific Safety Notes

Your weigh module meets the state of the art technology and complies with all recognized safety rules, however, certain hazards can arise.

Do not open the weigh module: It does not contain any parts which can be maintained, repaired or replaced by the user. If you ever have problems with your weigh module, contact your authorized METTLER TOLEDO dealer or service representative.

Observe instructions

Always operate and use your weigh module only in accordance with the instructions contained in the product documentation. The instructions for setting up your weigh module must be strictly observed.

If the weigh module is not used according to the product manuals, protection of the weigh module may be impaired and METTLER TOLEDO assumes no liability.

Staff safety

Use only METTLER TOLEDO accessories and peripheral devices, these items are designed to work optimally with your weigh module.

Explosion hazard

It is not permitted to use the weigh module in explosive atmospheres of gases, steam, fog, dust and flammable dust (hazardous environments).

3 WMF Weigh Modules

3.1 Specifications

Parameter		WMF204C	WMF303C
Maximum capacity	Nominal	220 g	320 g
Readability	Nominal	0.1 mg	1 mg
Measurement properties			
Specification temperatures		+10 ... +30 °C	
Specification humidity		20 ... 80 % rH	
Limit values			
Repeatability ¹⁾ (at nominal load)	Sd	0.13 mg	1 mg
Linearity deviation	Sd	0.4 mg	2 mg
Eccentricity deviation (test load)	Sd	1 mg (100 g)	2 mg (100 g)
Sensitivity offset ¹⁾ (test load)		0.8 mg (150 g)	2.0 mg (300 g)
Sensitivity temperature drift ³⁾		0.00015 % / °C x Rnt	0.00015 % / °C x Rnt
Sensitivity stability ⁴⁾		0.0025 % x Rnt	0.0025 % x Rnt
Typical values			
Repeatability	Typical	0.08 mg (200 g)	0.4 mg (300 g)
Eccentricity deviation (test load)	Typical	0.4 mg (200 g)	1 mg (300 g)
Sensitivity offset ¹⁾ (test load)	Typical	0.5 mg (150 g)	0.8 mg (300 g)
Minimum weight (according to USP)		160 mg	800 mg
Minimum weight (U = 1 %, 2 sd)		16 mg	80 mg
Dynamics			
Signal processing output rate	Max.	92/s	
Settling time ²⁾	Typical	< 0.4 s	

1) applies only after adjustment at nominal capacity with an OIML E2 weight

2) with a tolerance of 0.1 % of the applied load

3) weighing tests (loading and unloading) according to OIML R76 A.5.3


4) span stability test according to OIML R76 B4, applies only after internal adjustment every 8 hours


Rnt net weight (of sample)

sd standard deviation

3.2 Overview

WMF weigh module		
	1	Housing, stainless steel (1.4404-316L)
	2	Weighing platform, stainless steel (1.4404-316L)
	3	Status LEDs / Reed switch
	4	Base plate, stainless steel (1.4404-316L)

Connections at the bottom		
	5	M12-4, D coded
	6	Air inlet for cooling function
	7	Air outlet for cooling function
	8	Air connection for wash-down function
	9	Deaeration plug for wash-down function
	10	3 x M5 threaded holes, 6 mm deep, to mount the weigh module on a support platform
	11	Base plate sealing to seal the bottom plate on a support platform

Status LEDs		
	Module status	
	12	Status, bicolor green/red
	Industrial Ethernet status	
	13	NS / BF, bicolor green/red
	14	MS / SF, bicolor green/red
	Reed switch	
	15	Reed switch function (between the LED on the top and the one in the middle)

(13), (14): Refer to the Installation Manual for more details on the functions and meaning of the LEDs.

3.3 Scope of Delivery

All models are supplied by default with the following items:

- WMF weigh module incl. adapter weighing platform
- User manual (this document)
- Production certificate
- Declaration of Conformity
- Gasket bottom WMF

4 Operation

4.1 Electrical Connections

WMF weigh modules have an EtherNet/IP or PROFINET IO RT interface for the communication with the control systems.

Interface descriptions

EtherNet/IP	Adapter
PROFINET IO RT	Device, RT_CLASS_1, Conformance Class CC-B

Power supply of the weigh module

The weigh module can be energized with the Power over Ethernet (PoE).

- Power over Ethernet (PoE)
- Mode A (end span)
- Class 1 PD (below 3.84 Watt)
- According to IEEE Std. 802.3af

4.2 Module Status LED

The status of the weigh module is displayed with the bicolor LED on the very top (see pos. 12 in [Overview ▶ Page 6]).

Status	Meaning
Green	Normal operation
Red, blinking	Warning, see Installation Manual for further instructions
Red, solid	Error, see Installation Manual for further instructions

4.3 Environmental Conditions

WMF weigh modules can be operated within the following environmental conditions:

Temperature Range	Operating / Compensated	+10 ... +30 °C
	Allowable ambient	+5 ... +40 °C
Relative air humidity	Max. 80 % at 31 °C, decreasing linearly to 50 % at 40 °C, non-condensing	
Height above mean sea level	Max. 4'000 m (13'330 ft)	
Warm-up time	At least 45 minutes after power up	

NOTICE

Multi-line arrangement

Multi-line arrangement can only be implemented with the cooling function activated. For more details refer to the Installation Manual.

4.4 Ingress Protection

IP rating of the WMF weigh modules:

- IP44 in weighing configuration
- IP65 if the wash-down function is activated

4.5 Commissioning

Follow these steps to put the WMF weigh module into operation:

NOTICE

Damage to the weigh module

- 1 Always handle the weigh module with care.
- 2 Do not let the weigh module fall to the ground.

- 1 Install the weigh module in the place of use. For details refer to the Installation Manual.
- 2 Add custom-specific parts.
- 3 Switch on the power.
- 4 Observe the warm-up time (at least 45 minutes after power up).
- 5 Perform an adjustment before starting the weighing operation.

5 Maintenance

5.1 Cleaning

Dry Cleaning

- Use a damp cloth to clean the housing of the weigh module.
- Keep the area between the weighing platform and the upper part of the housing clean to ensure perfect operation of the weigh module.

Water Spray

NOTICE

Damage to the weigh module

- Make sure that the wash-down protection mechanism is activated when cleaning the weigh module by water spray.

Chemical Cleaning agents

WMF weigh modules are designed to resist common chemical cleaning procedures.

- Check the exact chemical resistance of the weigh module against the used chemical agent before starting with the cleaning process.

Important Steps after Cleaning

Before starting with the weight measurements:

- Wait until the weigh module cools down back to the operating temperature range.
- When cooled down, dry clean the surface.
- Perform an internal calibration and check the weighing performance of the weigh module.
- Adjust the weigh module, if necessary.

5.2 Calibration and Adjustment

As your weigh module is a precision measuring instrument, periodic maintenance is a prerequisite for perfect operation. The maintenance intervals depend on usage, ambient and environmental conditions. Maintenance work may only be performed by a METTLER TOLEDO service technician.

Checking the weighing performance

The precision of the weigh module is normally monitored by the tests function.

It is recommended, that the linearity, repeatability and other key figures of the weigh module are checked by a METTLER TOLEDO service technician.

Contact us about a service agreement tailored to your needs and budget.

5.3 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

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© Mettler-Toledo GmbH 02/2018
Mettler-Toledo GmbH